




CT2784WSA Laboratory Single Pure Water Faucet, Tin-Lined Brass, Panel Mounted, Self-Closing Control, 6" Rigid/Swing Vacuum Breaker Gooseneck



NOTES:

1. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.
2. FIXTURE HAS INTERIOR LINING OF PURE TIN.
3. FIXTURE HAS COLORTECH EPOXY POWDER COATED FINISH. SPECIFY COLOR (WHITE, GRAY OR TAN) WHEN ORDERING.

Measurements may vary ± 1/4".