



WaterSaver

ColorTech BT

BT2860I-231WSA

LABORATORY PUSH/TURN GAS VALVE FITTING, DECK MOUNTED SINGLE, POP-UP INDICATO DISCS

POWDER-COATED METAL HOODED
HANDLE WITH POP-UP INDICATORS
COLOR-CODED PER EN STANDARD

FORGED BRASS
VALVE BODY

REMOVABLE SEVEN
SERRATION HOSE END

G3/8 OUTLET

POSITIONING PINS
(SEE DETAIL A)

G1/2 MOUNTING SHANK
WITH LOCKNUT AND WASHER

5-11/16"
(145mm)

2-7/16"
(62mm)

Ø1-5/8"
(41mm)


2-3/8"
(60mm)

1-5/16"
(33.3mm)

Ø1/8" (3.2mm)
POSITIONING PINS
3/32" (2.4mm) EXPOSED LENGTH

DETAIL A

REPLACEABLE PUSH/TURN
CERAMIC DISC CARTRIDGE



MEASUREMENTS MAY VARY $\pm 1/4"$ (6mm).

NOTES:

1. FITTING IS FACTORY ASSEMBLED AND TESTED AT 250 PSI (17 BAR). MAXIMUM WORKING PRESSURE IS 150 PSI (10 BAR).
2. SPECIFY IF CLEANING FOR OXYGEN AND HIGH PURITY GASES IS REQUIRED.
3. FITTING HAS COLORTECH EPOXY POWDER COATED FINISH. SPECIFY COLOR (WHITE, GRAY OR TAN) WHEN ORDERING.