




# WaterSaver


## ColorTech BT

### BT2870-231WSA


LABORATORY FINE CONTROL NEEDLE VALVE FITTING, DECK MOUNTED SINGLE



REPLACEABLE  
STAINLESS STEEL SEAT



RENEWABLE FINE METERING  
STAINLESS STEEL FLOATING  
CONE



DETAIL A

#### NOTES:

1. FITTING IS CLEANED, DEGREASED AND PACKED FOR HIGH PURITY GAS SERVICE.
2. FITTING IS FACTORY ASSEMBLED AND TESTED AT 375 PSI (26 BAR). MAXIMUM WORKING PRESSURE IS 250 PSI (17 BAR).
3. FITTING HAS COLORTECH EPOXY POWDER COATED FINISH. SPECIFY COLOR (WHITE, GRAY OR TAN) WHEN ORDERING.

MEASUREMENTS MAY VARY  $\pm 1/4$ " (6mm).

WaterSaver Faucet Co. 312 666 5500 t  
701 West Erie Street 312 666 5501 f  
Chicago, IL 60654 wsflab.com

Drawing Number:

Revised: 022013-TEC