




BT3173-771WSA

PRESSURE REGULATOR FITTING, HIGH PURITY GASES, OVERHEAD CARRIER MOUNTED SINGLE

THIS FIXTURE IS FOR USE ON HIGH PURITY GAS SYSTEMS, INCLUDING OXYGEN. THE FOLLOWING GASES MAY BE USED WITH THIS FIXTURE: AIR, ARGON, CARBON DIOXIDE, HELIUM, NITROGEN, NITROUS OXIDE AND OXYGEN. DO NOT USE WHERE PRESSURE OR TEMPERATURE CAN EXCEED RATED OPERATING CONDITIONS.



MEASUREMENTS MAY VARY ±1/4" (6mm).