





BT7834

PURE WATER TAP, POLYPROPYLENE-LINED, PANEL MOUNTED, MANUAL CONTROL

BT7834SC

SAME AS ABOVE EXCEPT WITH SELF-CLOSING CONTROL



DETAIL A

NOTES:

1. MAXIMUM WORKING PRESSURE IS 100 PSI (7 BAR).
2. ALL COMPONENTS IN CONTACT WITH PURE WATER ARE INERT POLYPROPYLENE.
3. FITTING HAS COLORTECH EPOXY POWDER COATED FINISH. SPECIFY COLOR (WHITE, GRAY OR TAN) WHEN ORDERING.

MEASUREMENTS MAY VARY ±1/4" (6mm).

WaterSaver Faucet Co. 312 666 5500 t
 701 West Erie Street 312 666 5501 f
 Chicago, IL 60654 wsflab.com

Drawing Number:

Revised: 022113-TEC