



Stainless steel braided hose with
corrugated stainless steel inner core

Corrugated Stainless Steel Hose with Stainless Steel Overbraid

WaterSaver stainless steel/stainless steel hose assemblies are designed for use with our keyed quick connect fittings. Features of WaterSaver hose assemblies include the following:

- Hose has an inner core of corrugated stainless steel and a braided stainless steel exterior layer. Hose is thus resistant to cutting or abrasion. Hose is rated to 1000 PSI burst pressure, 300 PSI maximum working pressure.
- Braided stainless steel hose assemblies may be specially cleaned and packaged for high purity gas applications.
- End fittings are permanently attached to hose with a barb fitting and crimped ferrule. End connections swivel to prevent kinking or twisting of hose. Hoses are available with a variety of end connections for maximum flexibility.

Stainless Steel Hose Assemblies - Specify overall length of hose in inches (XXX) when ordering	
Item Number	End Connections
HSSS-XXX-3M	3/8" NPT MALE END CONNECTIONS
HSSS-XXX-3F	3/8" NPT FEMALE END CONNECTIONS
HSSS-XXX-3FM	3/8" NPT MALE X 3/8" NPT FEMALE

