



WaterSaver

INNOVATIVE PRODUCTS FOR RESEARCH


WaterSaver Faucet
701 W Erie St
Chicago, IL 60654

312 666 5500 TELEPHONE
312 666 5501 FACSIMILE
wsflab.com

L774VR

PRESSURE REGULATOR ASSEMBLY, HIGH PURITY GASES (2-125 PSIG), FUME HOOD MOUNTED,
VERTICAL INSTALLATION (RIGHT HAND)

THIS FIXTURE IS FOR USE ON HIGH PURITY GAS SYSTEMS, INCLUDING OXYGEN. THE FOLLOWING GASES MAY BE USED WITH THIS FIXTURE: AIR, ARGON, CARBON DIOXIDE, HELIUM, NITROGEN, NITROUS OXIDE AND OXYGEN. DO NOT USE WHERE PRESSURE OR TEMPERATURE CAN EXCEED RATED OPERATING CONDITIONS.



NOTES:

1. FIXTURE IS CLEANED, DEGREASED AND PACKED FOR HIGH PURITY GAS SERVICE.
2. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.
3. FIXTURE HAS A POLISHED CHROME PLATED FINISH.

MEASUREMENTS MAY VARY ±1/4" (6mm).

Drawing Number: _____

Revision Number: 020714-TEC