




WaterSaver


ColorTech

CT66VB-WSA

COMBINATION FIXTURE FOR COLD WATER AND GAS, 6" RIGID VACUUM BREAKER GOOSENECK



- NOTES:
1. WATER VALVE IS CERTIFIED BY CSA INTERNATIONAL TO COMPLY WITH ASME A112.18.1M AND CSA B125.1.
 2. VACUUM BREAKER IS CERTIFIED TO COMPLY WITH ASSE 1001.
 3. LABORATORY BALL VALVES ARE CERTIFIED BY CSA INTERNATIONAL TO COMPLY WITH ANSI Z21.15 AND CGA 9.1.
 4. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.
 5. FIXTURE MAY BE USED FOR COLD WATER AND NATURAL GAS, AIR OR VACUUM. UNLESS OTHERWISE SPECIFIED, FIXTURE IS FURNISHED WITH INDEX DISCS FOR GAS.
 6. FIXTURE HAS COLORTECH EPOXY POWDER COATED FINISH. SPECIFY COLOR (WHITE, GRAY OR TAN) WHEN ORDERING.



WATER SAVER RENEWABLE UNIT WITH REPLACEABLE STAINLESS STEEL SEAT

MEASUREMENTS MAY VARY ±1/4" (6mm).

WaterSaver Faucet Co. 312 666 5500 t
701 West Erie Street 312 666 5501 f
Chicago, IL 60654 wsflab.com

Drawing Number:

Revised: 012313-TEC