





CT681 Pure Water Faucet, Tin-Lined Brass, Deck Mounted, Self-Closing Control, 6" Rigid Gooseneck



TIN-LINED SELF-CLOSING WATER SAVER RENEWABLE UNIT WITH REPLACEABLE STAINLESS STEEL SEAT

NOTES:

1. FIXTURE HAS INTERIOR LINING OF PURE TIN.
2. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.
3. FIXTURE HAS COLORTECH EPOXY POWDER COATED FINISH. SPECIFY COLOR (WHITE, GRAY OR TAN) WHEN ORDERING.

Measurements may vary ± 1/4".