





WaterSaver

ColorTech BT

BT1714BH-WSA

LABORATORY MIXER TAP, PANEL MOUNTED, 6" RIGID/SWING GOOSENECK, WRIST BLADE HANDLES



WATER SAVER RENEWABLE UNIT WITH REPLACEABLE STAINLESS STEEL SEAT

MEASUREMENTS MAY VARY ±1/4" (6mm).

NOTES:

1. FITTING IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.
2. FITTING HAS COLORTECH EPOXY POWDER COATED FINISH. SPECIFY COLOR (WHITE, GRAY OR TAN) WHEN ORDERING.
3. IF CERAMIC DISC UNITS ARE DESIRED, ORDER BTC1714BH-WSA