



### BT3173-772WSA

PRESSURE REGULATOR FITTING, HIGH PURITY GASES (2-125 PSIG), OVERHEAD CARRIER MOUNTED SINGLE

THIS FIXTURE IS FOR USE ON ULTRA HIGH PURITY GAS SYSTEMS. THE FOLLOWING GASES MAY BE USED WITH THIS FIXTURE: AIR, ARGON, CARBON DIOXIDE, HELIUM, HYDROGEN, NITROGEN, NITROUS OXIDE AND OXYGEN. DO NOT USE WHERE PRESSURE OR TEMPERATURE CAN EXCEED RATED OPERATING CONDITIONS.

