





WaterSaver

ColorTech BT

BT5100F-226WSA

LABORATORY WATER VALVE FITTING, ANGLE PATTERN, PANEL MOUNTED



WATER SAVER RENEWABLE
UNIT WITH REPLACEABLE
STAINLESS STEEL SEAT

MEASUREMENTS MAY VARY $\pm 1/4"$ (6mm).

NOTES:

1. FURNISHED WITH THREADED ESCUTCHEON, MOUNTING SHANK, LOCKNUT AND WASHER FOR SECURING TO PANEL.
2. FITTING IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.
3. IF CERAMIC DISC UNIT IS DESIRED, ORDER BTC5100F-226WSA.
4. FITTING HAS COLORTECH EPOXY POWDER COATED FINISH. SPECIFY COLOR (WHITE, GRAY OR TAN) WHEN ORDERING.

WaterSaver Faucet Co. 312 666 5500 t
701 West Erie Street 312 666 5501 f
Chicago, IL 60654 wsflab.com

Drawing Number:

Revised: 021913-TEC