





WaterSaver

ColorTech BT

BT611BH

LABORATORY SINGLE TAP, DECK MOUNTED, 6" RIGID GOOSENECK, WRIST BLADE HANDLE



WATER SAVER RENEWABLE UNIT WITH REPLACEABLE STAINLESS STEEL SEAT

MEASUREMENTS MAY VARY ±1/4" (6mm).

NOTES:

1. FITTING IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.
2. FITTING HAS COLORTECH EPOXY POWDER COATED FINISH. SPECIFY COLOR (WHITE, GRAY OR TAN) WHEN ORDERING.
3. IF CERAMIC DISC UNIT IS DESIRED, ORDER BTC611BH.

WaterSaver Faucet Co. 312 666 5500 t
 701 West Erie Street 312 666 5501 f
 Chicago, IL 60654 wsflab.com

Drawing Number:

Revised: 022013-KJS