



L2794WSA Laboratory Single Pure Water Faucet, Tin-Lined Brass, Panel Mounted, 6" Rigid/Swing Gooseneck



TIN-LINED WATER SAVER RENEWABLE UNIT WITH REPLACEABLE STAINLESS STEEL SEAT

NOTES:

1. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.
2. FIXTURE HAS INTERIOR LINING OF PURE TIN.
3. FIXTURE HAS A POLISHED CHROME PLATED FINISH.

Measurements may vary $\pm 1/4"$.