


- ☐ **BT081** Bench Mounted Gooseneck Fitting, Rigid Gooseneck
- ☐ **BT084** Bench Mounted Gooseneck Fitting, Rigid/Swing Gooseneck



**Application:** Bench or countertop mounted gooseneck fitting for water.

**Body:** Brass bar stock.

**Gooseneck:** 150mm (6") brass. BT081 has rigid gooseneck. BT084 has convertible rigid/swing gooseneck.

**Escutcheon:** Molded nylon three-piece escutcheon color coded and marked per EN 13792.

**Inlet:** Furnished with G1/2 mounting shank, locknut and washer. Available with choice of inlet connections (see below).

**Outlet:** G3/8 female outlet with removable anti-splash serrated hose end.

**Quality Assurance:** Fitting is fully assembled and factory tested prior to shipment. 10 bar (150 PSI) test pressure and maximum working pressure.

Ordering Information

	Inlet Connection Options			
	G1/2 Male	G1/4 Female	1/4 NPT Female	10mm Push/Fit Tube
Rigid Gooseneck	<input type="checkbox"/> BT081	<input type="checkbox"/> BT081-G2F	<input type="checkbox"/> BT081-2F	<input type="checkbox"/> BT081-M10PF
Rigid/Swing Gooseneck	<input type="checkbox"/> BT084	<input type="checkbox"/> BT084-G2F	<input type="checkbox"/> BT084-2F	<input type="checkbox"/> BT084-M10PF

Product Finish			
<input type="checkbox"/> WHT	Gloss white powder coated finish	<input type="checkbox"/> PCL	Polished chrome plated finish with clear epoxy coating
<input type="checkbox"/> GRY	Gloss light gray powder coated finish	<input type="checkbox"/> SCC	Satin chrome plated finish with clear epoxy coating
<input type="checkbox"/> GRD	Gloss dark gray powder coated finish	<input type="checkbox"/> SNC	Satin nickel plated finish with clear epoxy coating
<input type="checkbox"/> STM	Starburst metallic powder coated finish		

**Note:** Refer to "Product Finishes" catalog for complete information on finishes.

Media

☐