



### ○ L2880 Needle Valve, Straight Pattern



REPLACEABLE  
STAINLESS STEEL SEAT



STAINLESS STEEL  
REPLACEABLE FLOATING NEEDLE

#### NOTES:



1. VALVE IS CERTIFIED BY CSA INTERNATIONAL TO COMPLY WITH ANSI Z21.15 AND CGA 9.1.
2. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED AT 250 PSI NITROGEN PRESSURE. MAXIMUM WORKING PRESSURE IS 150 PSI.
3. SPECIFY IF CLEANING FOR OXYGEN AND HIGH PURITY GASES IS REQUIRED.
4. VALVE HAS A POLISHED CHROME PLATED FINISH.

Measurements may vary  $\pm 1/4"$ .