





WaterSaver

VR5600WSA

VANDAL-RESISTANT COMBINATION FIXTURE, HOT/COLD WATER AND GAS



WATER SAVER RENEWABLE UNIT WITH REPLACEABLE STAINLESS STEEL SEAT

MEASUREMENTS MAY VARY $\pm 1/4$ " (6mm).

WaterSaver Faucet Co. 312 666 5500 t
701 West Erie Street 312 666 5501 f
Chicago, IL 60654 wsflab.com

Drawing Number:

Revised: 110413-TEC



NOTES:

1. WATER VALVE IS CERTIFIED BY CSA INTERNATIONAL TO COMPLY WITH ASME A112.18.1M AND CSA B125.1.
2. LABORATORY BALL VALVES ARE CERTIFIED BY CSA INTERNATIONAL TO COMPLY WITH ANSI Z21.15 & CGA 9.1.
3. ALL THREADED CONNECTIONS ARE SECURED WITH SET SCREWS OR CEMENT.
4. LABORATORY BALL VALVES HAVE INTEGRAL CHECK VALVES TO PREVENT BACKFLOW.
5. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.
6. FIXTURE HAS A POLISHED CHROME PLATED FINISH.