




WaterSaver

L3173-366-758WSA

PRESSURE REGULATOR FIXTURE, HIGH PURITY GASES (2-125 PSIG), PANEL MOUNTED SINGLE, FLOWSTOP, MOUNTING SHANK

THIS FIXTURE IS FOR USE ON HIGH PURITY GAS SYSTEMS, INCLUDING OXYGEN. THE FOLLOWING GASES MAY BE USED WITH THIS FIXTURE: AIR, ARGON, CARBON DIOXIDE, HELIUM, NITROGEN, NITROUS OXIDE AND OXYGEN. DO NOT USE WHERE PRESSURE OR TEMPERATURE CAN EXCEED RATED OPERATING CONDITIONS.



NOTES:

1. FIXTURE IS CLEANED, DEGREASED AND PACKED FOR HIGH PURITY GAS SERVICE.
2. FIXTURE IS FULLY ASSEMBLED AND FACTORY TESTED PRIOR TO SHIPMENT.
3. FIXTURE HAS A POLISHED CHROME PLATED FINISH.

MEASUREMENTS MAY VARY $\pm 1/4"$ (6mm).